

F16

Easy to attach / detach powerchair add on kit...



Easy to attach / detach powerchair add on kit...

One hand single motion detachment

Includes easy disassembly of battery packs for ease of transportation.

Easy fitting to a wide range of manual chairs

Easily fitted to a broad range of frame sizes; 34 to 48cm (13.5 to 19").

Performance

Stable rear wheel adjustment for better stability and handling.

Powerful 25Ah gel batteries capable of up to 9.3 miles (15km) at 3.7mph (6kph).

- Joystick controlled operation for those with restricted hand mobility
- Compact for use indoors, with a tight turning circle

Fits to:

- Neon, Neon Swing-Away, Youngster 3, RXS, GPV, RubiX, and Easymax manual chairs



One hand fixation and detachment



25Ah gel batteries



Removes for transport / storage in seconds

Technical Specifications

Total width F16:	51 - 65cm (20 - 25.5")	Battery:	2 x 25Ah
Max safe slope:	10 degrees	Max load:	125kg (19.6 stone)
Max. range:	15km (9 miles)	Weight (Power F16 excl. batteries):	23kg (50.7lbs)
Speed:	3.7mph (6kph)	Weight (Power F16 incl. 2 batteries):	43.8kg (96.5lbs)

Find out more today!

Visit www.quickiepower.com for more information, to find your local dealer or to request a demonstration.

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Rumba and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairs



Sunrise Medical Limited
Thorns Road
Brierley Hill
West Midlands
DY5 2LD
Tel: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89
www.SunriseMedical.com