# Invacare Action 4<sub>NG</sub>







## Where stability meets style!

Invacare  $Action^4_{NG}$  offers a flexible range of models to accommodate different customer needs and is particularly suitable for taller or heavier people.

With a frame designed for improved stability and comfort, the stylish  $\mathbf{Action^a}_{NG}$  offers the same cross compatibility of parts with other models in the  $\mathbf{Action^a}_{NG}$  family allowing customers to easily maintain and update their wheelchairs when required.





#### Easy rolling and maneuverability

The frame design of **Action**<sup>4</sup><sub>NG</sub> helps improve stability resulting in easier wheeling and manoeuvring by maximising the customer's position in the chair.

It is easy to adapt the **Action**<sup>4</sup><sub>NG</sub> to other types of wheeling characteristics to suit various needs by changing the rear wheel position and altering the seat position. It features optional seat adjustments to aid postural support for day long use.

# Compatibility with other Invacare wheelchairs

Almost all parts of the **Action**<sup>4</sup><sub>NG</sub> frame are interchangeable with other models of the **Action** family and even with other Invacare wheelchairs, including *Rea*<sup>\*\*</sup>. All parts such as armrests, legrests, wheels etc, are compatible with each other. The **Action**<sup>4</sup><sub>NG</sub> is also compatible with most of the **Alber** product range - please see the **Alber** brochures for further details.

#### Easy transport with low weight

By removing the rear wheels and folding the frame, the  $Action^4_{NG}$  easily slides in behind the front seat or slips into the boot of a car. Without rear wheels, the  $Action^4_{NG}$  weighs 8.3 kg. Simple handling and folding makes for great ease of transport.

Even the most basic  $Action^4_{NG}$  model covers customer needs, as almost everything, from the height of

the backrest, armrest and footrest is adjustable.

From the  $Action_{NG}^4$  Tall (125 kg) Transit model for use with an attendant, to the self-propelling model or the groundbreaking  $Action_{NG}^4$  Heavy Duty version (with a user weight of 160 kg) the  $Action_{NG}^4$  offers total flexibility in a chair.

#### Features and Options



Available as either a standard Transit version or Transit HD.



## Adjustable backrest upholstery

As an option, the backrest can be equipped with tensionadjustable straps, shaping the backrest to individual needs.



**Seat-depth adjustability**The seat depth has up to 50 r

The seat depth has up to 50 mm adjustment with the upholstery.



**Backrest bar** 

With tension adjustable upholstery a backrest bar is available for extra rigidity.

#### **Features and Options**



#### Hangers and legrests

Hangers come as angle-adjustable, fixed 80 degree angle and amputee version. Use depth-adjustable, fixed or one piece footrest with your hanger.



#### Optional rear wheel positions

New active position plate offers better balancing possibilities for superior driving performance.



#### Headrest

A new range of head supports are available offering a greater range of adjustment and positioning.



#### **Hub brakes**

Hub brakes available as an option, helping assistants to stop the wheelchair and reduce speed.



#### Anti-tip device

Prevents wheelchair from tipping over and is height adjustable to suit the different seat heights available.



#### **Amputee support**

The amputee support is both height and angle adjustable.



#### **Heavy Duty Version 160 kg**

A wider choice of seat widths (505, 555 & 605 mm), reinforced chassis and a greater carrying capacity (maximum user weight 160 kg).



### Hemi version

Dual handrim propulsion and hemi armrest are available.



400 / 450 / 500 mm



410 / 435 / 460 / 380 - 460 / 390 - 510 mm\* 485 / 510 mm



200 - 230 - 260 / 190 - 295 mm



Invacare Action4<sub>NG</sub> Invacare Action HD

455 / 480 mm

505 / 555 / 605 mm 400 / 450 / 500 mm

410 / 435 / 460 / 485 / 510 mm

380 - 460 / 390 - 510 mm\*

230 - 260 - 280 mm 350 - 480 mm 200 -190 - 295 mm



0° - 80°

Seat width



900 - 1030 mm



1050 - 1110 mm



830 - 900 mm

Invacare Action⁴<sub>NG</sub>

Invacare Action⁴<sub>NG</sub> HD

0° - 30°\* 7° / 0° - 30°\*

+ 200 mm Seat width + 200 mm

900 - 1030 mm

1050 - 1110 mm 830 - 900 mm



305 mm



15.2 kg



125 kg



8.3 kg

8.3 kg



Crash tested ISO 7176-19\*\*

Crash tested ISO 7176-19\*\*



Invacare Action4<sub>NG</sub>

Invacare Action HD

305 mm

17 kg

160 kg\*\*

\* Reclining backrest

## Frame colours



Dark Red



Slate Grey



Ocean Blue



Pearl White



Anthracite Black



Black Nylon



Please notice colors may vary slightly from those displayed above.

Invacare Limited Pencoed Technology Park - Pencoed -Bridgend CF35 5AQ - United Kingdom Tel: +44 1 656 776 222 Fax: +44 1656 776 220

E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com www.invacare.co.uk

Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS - County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie

© 2013 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Action4NG - UK - 03/2013 1528015



<sup>\*\*</sup> Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.