

Invacare® | Storm⁴ Max



Much more than just a larger power chair!

The new *Invacare Storm⁴ Max* is designed specifically for the needs of clients with a larger bodyshape. Its longer chassis supports the shifted centre of bodymass for improved weight distribution and swift mobility. Extra strong arm- and legrest support the clients needs to push themselves back to re-position or do some type of pressure relief.

The **G-Trac[®]** option is also available, realizing the ultimate driving performance and traction control.

The stylish **Storm⁴ Max** completes the **Storm⁴** family and conserves the trouble-free servicing.



Yes, you can.®



Features and Options

Tilt and recline

For individual comfort and positioning the **Storm⁴ Max** offers powered tilt and recline options. The backrest recliner has a range of 0° - 25° with a high pivot point in order to reduce shearing forces.

Ease of use

Just as the entire *Invacare Storm⁴* family is characterised by its ease of use, so is the **Max**. Batteries and electronics are easily accessible, making servicing simple. The depth of the seat can be simply adjusted in 3 steps which are indicated on the seat frame. An adapter to further reduce the seat depth is offered as an option.



Armrest

The strong and stylish armrest can be mounted in several positions as well as adjusted in width.



Footboard

Flip-up footboard supports various positions of the feet. It's strong enough to allow the clients to push themselves back to reposition or do some type of pressure relief.



Calf pads

The extra wide Max calf pads are individually adjustable in depth and angle to give optimum support.

Technical data

Storm⁴ Max	490 - 590 mm 550 - 650 mm 610 - 710 mm	560, 595, 630 mm ¹	480, 510 mm	570 mm ^{2,3} , 630 mm ³	360 - 480 mm	0°- 16° ⁴ , 0°- 9° ⁵	0 - 25°	1390 mm
Storm⁴ Max	2080 mm	60 mm	6°(10,5%)	350W	31 km ⁶	6 km/h	220 kg	

Standard colours



Blue (RAL 5026) Black (RAL 9005) Champagne colour (0314)

Colour options



Transparent red (0200) White (RAL 9016) Olive green (RAL 6003)

- 1/ Adapter to reduce seat depth by 130 mm available
- 2/ For small seat frame (490 mm) only
- 3/ Additional 50 mm of height adjustment
- 4/ Powered seat tilt
- 5/ Man. settable seat tilt
- 6/ according to ISO 7176-4:2008



More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.